



# Gas Exchange Instrument theory

How the LI-6800 works and  
what it can tell you


Copyright © 2019 LI-COR, Inc. All rights reserved.  
Contact LI-COR for permission to redistribute.

Used by permission from LICOR (November, 2023)

**LI-COR**<sup>®</sup>

# How do we measure photosynthesis?


- Photosynthesis
- Gas exchange
- Assimilation



# How do we measure photochemistry?

Gas exchange

Chlorophyll  
Fluorescence



# How do we measure gas exchange?

- First we need to measure the concentration of one or more gases

$$A'_{(\nu)} = 1 - e^{-k_{(\nu)}wl}$$

$w$  = density of the absorber

$l$  = path length

$\nu$  = wavelength

$A'$  = light absorption

$k(\nu)$  = absorption coefficient

# How do we measure gas exchange?

- Infrared Absorption by Gases

$$\alpha + T = 1$$

$$\alpha = 1 - T$$


$$\alpha = 1 - \Phi_A / \Phi_{NA}$$

$\alpha$  = absorptance


$T$  = transmittance

$\Phi_A$  = transmittance in  
absorbing region

$\Phi_{NA}$  = transmittance in  
non-absorbing region



# Basic optical differential IRGA




$$\alpha = 1 - \Phi_A / \Phi_{NA}$$

$\alpha$  = absorptance  
 $\Phi_A$  = transmittance in absorbing band  
 $\Phi_{NA}$  = transmittance in non-absorbing band

# Gas exchange systems

- Closed-transient
- Compensating



# Gas exchange systems


- Closed-transient
- Compensating
- Open flow through






# System flow path


Traditional open system



LI-6400/XT




LI-6800




# Basic parameters *computed*

- Fluxes from mass balance
  - A – Assimilation rate
  - E – Apparent Transpiration rate
- Additional parameters
  - $g_{sw}$  – Stomatal conductance to water vapor
  - $C_i$  – Inter-cellular CO<sub>2</sub> concentration

# Mass balance in an open system



# Mass balance for gas exchange



$$sf_x = u_{out}x_{out} - u_{in}x_{in} \quad \text{Steady-state gas exchange model}$$

$s$  = Leaf surface area

$f_x$  = Flux of gas  $x$


$u_{in}$  = Molar air flow entering the chamber

$u_{out}$  = Molar air flow exiting the chamber

$x_{in}$  = Mole fraction of  $x$  entering chamber

$x_{out}$  = Mole fraction of  $x$  exiting the chamber

# Mass balance for gas exchange



$$sf_x = u_{out}x_{out} - u_{in}x_{in}$$

Steady-state gas model


$$sf_x = u_{out}x_{out} - u_{in}x_{in} + \rho v \frac{dx}{dt}$$

Dynamic model

When steady-state assumptions are met  $\rho v \frac{dx}{dt}$  evaluates to 0

# Mass balance for gas exchange

Substituting in for measured values and rearranging the dynamic model...



Flow out of the chamber is not measured. It is accounted for by assuming  $u_{out} = u_{in} + sE$ , and  $E$  is the only significant flux.  $E$  is generally  $10^3$  or  $10^4 > A$ .

$$E = \frac{u_{in}}{s(1 - w_{out})} \left( w_{out} - w_{in} + \frac{\rho v}{u_{in}} \frac{dw_{out}}{dt} \right) \quad A = \frac{u_{in}(1 - w_{in})}{s} \left( c'_{in} - c'_{out} + \frac{\rho v}{u_{in}} \frac{dc'_{out}}{dt} \right)$$

$E$  = Apparent transpiration


$A$  = Carbon assimilation

$w_x$  = Water vapor mole fraction

$c'_x$  =  $CO_2$  dry mixing ratio =  $\frac{c_x}{1 - w_x}$

# Mass balance for gas exchange

Substituting in for measured values and rearranging the dynamic model...



Flow out of the chamber is not measured. It is accounted for by assuming  $u_{out} = u_{in} + sE$ , and  $E$  is the only significant flux.  $E$  is generally  $10^3$  or  $10^4 > A$ .

$$E = \frac{u_{in}}{s(1 - w_{out})} \left( w_{out} - w_{in} + \frac{\rho v}{u_{in}} \frac{dw_{out}}{dt} \right)$$

$$A = \frac{u_{in}(1 - w_{in})}{s} \left( c'_{in} - c'_{out} + \frac{\rho v}{u_{in}} \frac{dc'_{out}}{dt} \right)$$

$E$  = Apparent transpiration

$A$  = Carbon assimilation


$w_x$  = Water vapor mole fraction

$c'_x$  =  $\text{CO}_2$  dry mixing ratio =  $\frac{c_x}{1 - w_x}$

At steady state, the transient terms are zero and the equations are equivalent to what has always been used in the LI-6400/XT and LI-6800

# Steady state versus dynamic

Example data from an empty chamber (Assimilation = 0)




Following a step change in CO<sub>2</sub> the steady-state model gives a false flux driven by washout.

The dynamic model accounts for this and still shows zero assimilation.



# Accounting for dilution



$$A \gg \frac{F(C_r - C_s)}{100S}$$

$$A = \frac{F(C_r - C_s) \frac{1000 - W_r}{1000 - W_s}}{100S}$$

# More on water corrections...



APPLICATION NOTE

## The Importance of Water Vapor Measurements and Corrections

Application Note #129

Water vapor is known to influence the measurement of carbon dioxide by infrared gas analysis in several ways, which can lead to significant measurement errors. Spectral cross-sensitivity due to absorption band broadening, and inherent instrument cross-sensitivity can both cause overestimations of CO<sub>2</sub> mole fraction in samples containing water vapor when their effects are not accounted for. Dilution of samples by the addition of water vapor may not be important when measuring actual CO<sub>2</sub> mole fractions, but can lead to significant errors in flux measurements. In this note we describe the basis of each of these three processes and discuss how each can affect the measurement of CO<sub>2</sub>.


### Infrared Absorption by Gases and Absorption Band Broadening

At the sub-molecular scale the positions of atoms within molecules are not entirely fixed. As they move they stretch and bend their bonds, creating vibrations within the molecule. Energy differences between the possible vibrational states that result from this movement make it possible for the molecule to absorb infrared radiation. Changes in angular momentum as the molecule rotates about its axis can also cause infrared radiation to be absorbed. The energy differences resulting from changes in the vibrational and rotational states of the molecule cause fluctuations in its dipole moment. These oscillations interact with the alternating electrical field of electromagnetic radiation and if the frequencies of oscillation match, the radiation will be absorbed by the molecule.

Since the frequency of oscillation for electromagnetic radiation ( $\nu$ ) is the inverse of its wavelength ( $\lambda=1/\nu$ ) and the various energy states that lead to its absorption by a molecule are the result of the molecule's structure, the absorption of infrared radiation is both wavelength and absorber species dependent. For a given molecular species, absorption of infrared radiation will occur in bands at various wavelengths across the infrared region of the spectrum (Figure 1A). Each of these absorption bands is comprised of individual absorption lines (Figure 1B) that result from rotational transitions, and have a generally Lorentzian line shape described by

$$k(\nu) = \frac{S}{\pi} \frac{\alpha}{(\nu - \nu_0)^2 + \alpha^2} \quad (1)$$

where  $k(\nu)$  is the absorption coefficient  $k$  of radiation at frequency  $\nu$ ,  $S$  is the line strength,  $(\nu - \nu_0)$  is the change in frequency across the absorption line and  $\alpha$  is the absorption line half width.



**Figure 1: Absorption spectrum of CO<sub>2</sub>.** Carbon dioxide absorption data from HITRAN96 (Rothman et al. 1998). A. The absorption spectrum of CO<sub>2</sub> in the infrared region of the electromagnetic spectrum. B. The absorption of CO<sub>2</sub> in the 4.2  $\mu$ m region.

For infrared gas analysis with LI-COR gas analyzers, we are interested in measuring the total absorption of infrared light across a given waveband to determine absorber species concentration. The total absorption measured by the analyzer,  $A$  is in principle the result of integrating the absorption across the entire band as shown by

$$A = \int A'(\nu) d\nu \quad (2)$$

where  $A'(\nu)$  is the absorption of infrared radiation at a given radiation frequency, and is given by Beer-Lambert's Law

$$A'(\nu) = 1 - \exp[-k(\nu)wl] \quad (3)$$


where  $w$  is the absorber concentration and  $l$  is the optical path length (Burch and Williams 1964). These relationships form the basic principles behind

- Explains basis of broadening and effective pressure
- Derive dilution corrections
- How IRGAs work


# Matching IRGAs

LI-6800

Normal operating mode:




Match mode:



# Basic parameters *computed*

- Fluxes from mass balance
  - A – Assimilation rate
  - E – Transpiration rate
- Additional parameters
  - $g_{sw}$  – Stomatal conductance to water vapor
  - $C_i$  – Intercellular CO<sub>2</sub> concentration

# What else can we determine with gas exchange?



# What else can we determine with gas exchange?

- Fick's First Law

$$J_j = -D_j \frac{\partial c_j}{\partial x} = g_j \Delta c_j$$

$J_j$  = flux

$D_j$  = diffusivity coefficient

$\partial c_j / \partial x$  = change in concentration

$g_j$  = conductance

$\Delta c_j$  = concentration gradient

# What else can we determine with gas exchange?

$$\cancel{E} \approx g_{total}^{H_2O} (\cancel{W_i} - \cancel{W_a})$$

$$\cancel{A} \approx g_{total}^{CO_2} (\cancel{C_a} - C_i)$$

- Measure  $E$  &  $w_a$  ( $w_a = W_s$ )
- Measure leaf temperature
- Calculate  $w_i$
- Solve for  $g_{total}^{H_2O}$
- Measure  $A$  &  $c_a$  ( $c_a = C_s$ )
- $g^{CO_2} = g^{H_2O}/1.6$
- Solve for  $c_i$

# What else can we determine with gas exchange?


## Ohm's Law Analogy

$$r_{total} = r_{bl} + \frac{1}{\frac{1}{r_s} + \frac{1}{r_c}} + r_{mes}$$

Assumptions:


- end point of diffusion path is mesophyll surface
- cuticular resistance is near infinite

$$r_{total} \gg r_{bl} + r_s$$





# Boundary layer




$$S^{bl} = 4 \sqrt{\frac{l}{v}}$$

$S^{bl}$  = boundary layer thickness



$v$  = air velocity

$l$  = leaf width

# Boundary layer conditions



# Boundary layer conditions



# What does this mean for making measurements?

- For H<sub>2</sub>O: Target ambient boundary layer conditions!

$$\text{RH (\%)} = \frac{e_c}{e_{T_{air}}} * 100$$

$$e_c \text{ (kPa)} = W_s P$$


$$\text{VPD}_{leaf} \text{ (kPa)} = e(T_{leaf}) - e_c$$

On the LI-6800 use constant VPD (1 to 1.5 kPa)

# What does this mean for interpreting the data?

- $E$  = Real count of  $H_2O$  molecules leaving the leaf

$$g_s = \frac{E}{(w_i - w_a)}$$



# What does this mean for interpreting the data?

- Water Use Efficiency
- Instantaneous versus Intrinsic

$$W_t = \frac{A}{E} = \frac{A}{g_s (w_i - w_a)} = \frac{A}{g_s D_a} = W_g D_a$$


$$W_g = \frac{A}{g_s}$$

Seibt et al. 2008. *Carbon isotopes and water use efficiency: sense and sensitivity*. *Oecologia*

# Measurements fall in two categories

- Response Curve – imposed conditions
  - Light Response Curve
  - CO<sub>2</sub> Response Curve
- Survey
  - *In situ* picture under ambient conditions

# When to make measurements?






# A-C<sub>i</sub> and A-Q curves

- Response to [CO<sub>2</sub>]
- Response to light

# Parameters from AQ Response

$$A_{\text{sat}} = 34.28 \pm 0.778 \quad r^2 = 0.9995$$


$$R_d = -1.492 \pm 0.117$$



- $A_{\text{sat}}$  - photosynthesis at saturating light
- $R_d$  - Dark respiration rate
- LCP - Light compensation point
- $\phi$  - Quantum yield

# Parameters from $A_c$ response

- $V_{c,max}$  – velocity of carboxylation
- $J_{max}$  – electron transport for RuBP regeneration
- $c_i$  inflection – transition from  $V_{c,max}$  to  $J_{max}$
- $I$  – stomatal limitation to photosynthesis
- $g_m$  – mesophyll conductance
- $\Gamma^*$  –  $CO_2$  compensation point
- TPU – trios-phosphate utilization




# Considerations for $A_c$ measurements

- Saturating light
  - $\text{CO}_2$  effect on photosynthesis will be confounded with the effects of light
- Environmental Control
  - **Reference** vs. Sample
  - Humidity will alter stomatal aperture
  - Temperature control for enzyme kinetics
- Match at every point

# Using short time steps is important

- Rubisco deactivation
- Stomatal closure



# AC<sub>i</sub> curve fitting tools

PCE\_1710\_ACi curve fitting model.xls

Home Layout Tables Charts SmartArt Formulas Data Review

Edit Font Arial 10 Alignment Number General Format Cells Themes

Instructions for use  
 1. Enter data in white cells, delete any extra "A" values  
 2. Estimate limiting factors (1= rubisco, 2= RuBP regeneration, 3= TPU)  
 (assign at least one point to limitation 3, enter 0 to exclude points)  
 3. Press the "Solve" button  
 4. Adjust limiting factor if needed (use 0 to disregard a data pair)  
 5. Press the "Save" button to save to your computer and or  
 6. Cut and paste outputs if desired

Please enter your values  
 T leaf 28 °C  
 P<sub>atm</sub> 101 kPa  
 O<sub>2</sub> 21 kPa

Make no changes here

Estimate Limiting	Enter A	Enter either		Calculated C <sub>c</sub>	Limitations			Error terms	
		c, ppm*	C, Pa		Ac	Aj	At	0.480	0.110
1	-3.27	20.7	2.09	2.40	-3.62	-7.70	33.30	0.125	
1	3.49	77.5	7.83	7.50	4.69	8.01	33.30	0.291	
1	10.8	135	13.64	12.81	10.87	16.07	33.30	0.005	
1	17.8	197	19.90	18.21	17.86	21.32	33.30	0.059	
0	23.5	266	26.87	24.64	24.36	25.14	33.30		
0	27.6	344	34.74	32.12	31.31	28.06	33.30		
2	30.5	428	43.23	40.33	37.95	30.23	33.30	0.073	
2	31.7	517	52.22	49.21	44.18	31.87	33.30	0.031	
2	33	606	61.21	58.07	49.62	33.08	33.30	0.006	
3	33.3	698	70.50	67.34	54.60	34.04	33.30	0.000	
3	33.3	791	79.89	76.73	59.06	34.79	33.30	0.000	

\* If you enter c<sub>i</sub>, C<sub>i</sub> will be calculated, do not delete equations in the C<sub>i</sub> column

Use solver to minimize this sum of squares: **0.590**

Outputs

	@ T leaf	@ 25 °C
V <sub>c</sub> max	130	100 μmol m <sup>-2</sup> s <sup>-1</sup>
J	167	140 μmol m <sup>-2</sup> s <sup>-1</sup>
TPU	11.4	9.7 μmol m <sup>-2</sup> s <sup>-1</sup>
R <sub>g</sub> *	0.93	0.77 μmol m <sup>-2</sup> s <sup>-1</sup>
g <sub>m</sub> *	10.53	8.64 μmol m <sup>-2</sup> s <sup>-1</sup> Pa <sup>-1</sup>

\* R<sub>g</sub> is constrained to be >0 and g<sub>m</sub> is constrained to be 30 or less

To cite this estimating utility:  
 Sharkey, T.D., Bernacchi, C.J., Farquhar, G.D.,  
 Singaas, E.L. (2007) In Practice: Fitting photosynthetic  
 carbon dioxide response curves for C3 leaves. *Plant, Cell  
 & Environment* 30:XXX-XX  
 Version 2007.1  
[Link to paper](#)

To refresh your copy and for updates to this application please visit [www.blackwellpublishing.com/plantsci/pcecalculation/](http://www.blackwellpublishing.com/plantsci/pcecalculation/)

A/C<sub>i</sub> curve



Sharkey et al. 2015. What gas exchange data can tell us about photosynthesis. *Plant, Cell and Environment*

# AC<sub>i</sub> curve fitting tools

- Leafweb.ornl.gov
  - AC<sub>i</sub> curve fitting and data sharing platform
  - Developed and hosted by Oak Ridge National Laboratory
- Plantecophys R package
  - Duursma RA. 2015. Plantecophys – *An R package for analysing and modeling leaf gas exchange data*. PLoS ONE

# Other chamber options for the LI-6800

Our discussion of  $g_s$  and  $c_i$  thus far has been based on assumptions that break down at the whole plant/canopy scale!





Questions?






# Selecting the right chamber?


- Chamber aperture - bigger can be better!
- Light intensity control
- What data do you need?




# LI-6800 Console flow path:

Rapid & stable software control of CO<sub>2</sub> and H<sub>2</sub>O







# Leaf photosynthesis and carbohydrate dynamics of soybeans grown throughout their life-cycle under Free-Air Carbon dioxide Enrichment

A. ROGERS<sup>1,4</sup>, D. J. ALLEN<sup>2,3,\*</sup>, P. A. DAVEY<sup>1,2,†</sup>, P. B. MORGAN<sup>2</sup>, E. A. AINSWORTH<sup>1</sup>, C. J. BERNACCHI<sup>2,3</sup>, G. CORNIC<sup>5</sup>, O. DERMODY<sup>2</sup>, F. G. DOHLEMAN<sup>1</sup>, E. A. HEATON<sup>1</sup>, J. MAHONEY<sup>4</sup>, X.-G. ZHU<sup>2</sup>, E. H. DELUCIA<sup>2</sup>, D. R. ORT<sup>1,2,3</sup> & S. P. LONG<sup>1,2</sup>

